



High capacity washer disinfectors
FOR CSSD



For more than 30 years Smeg Instruments represents a symbol of excellence "Made in Italy". It offers technologically advanced products to the professional washing and thermal disinfection world.

Fully customisable solutions, flexible and efficient services, safe and certifiable results thanks to an integrated structure from design to after-sales service.

Products enriched by know-how gained from working with professionals in more than thirty countries all over the world.





A GLOBAL BRAND

Smeg is an Italian family company based in Guastalla (RE), with 4 production plants in Italy, and has been manufacturing household appliances since 1948. Smeg operates worldwide through subsidiaries, representative offices and an extensive network of distributors. Vision, intuition, international scope: the Smeg Group, thanks to a corporate culture that places extreme attention to quality, technological content and attention to energy saving, is recognised internationally as one of the best *Made in Italy* companies.

EXCELLENCE MADE IN ITALY

Thanks to its long and consolidated experience in the in-house production of technologies for the treatment of medical instruments, Smeg is able to guarantee its customers perfect thermal disinfection quality. Washer disinfectors and accessories specially designed to wash simple surgical instruments up to the delicate specific instruments for the robotic surgery.

Washer disinfectors designed and produced by Smeg, in Italy, are done so with our customers specific requirements in mind: compact, fast, ergonomic and reliable.

COMPACT WITH HIGH CAPACITY



The most compact washing disinfectors, allowing high volume application in reduced space, increased production incorporated within an ergonomic handling of the load.

FLEXIBILITY



Perfect treatment of different instruments. One rack, many applications and less costs.

FAST CYCLES



Complete washing and disinfection cycles performed in 35'. Less time, money and energy, lower running costs and working hours.

CERTAINTY OF RESULTS

Quality

Smeg's Quality System complies with 2017/745 (MDR) and ISO 9001:2015 and ISO 13485:2016 standards, as attested by the relevant certificates.

All washer-disinfectors, classified as class IIb medical devices, are CE certified and comply with EN 15883.

A0 parameter

Smeg washer-disinfectors offer a wide choice of programmes with preset or selectable A0 values depending on the type of instruments to be treated.

The value of A0 displays in real time and/or is printed out.

Traceability

The Smeg washer-disinfector range offers all connection and data storage options for traceability in line with the international EN ISO 15883 standards.

quality

safety

saving

WASHER DISINFECTOR WD6010

Washer disinfector with single or double drop down door in glass. Excellent washing capacity: up to 12 DIN baskets. The drop down door allows the customer to load and unload in an ergonomic way, and with our unique inserts and rack design adding safe and easy loading without any heavy lifting. Also featuring our 35' TURBO thermal disinfection cycle.

12
DIN

SINGLE OR DOUBLE
DROP DOWN DOOR




INTEGRATED PRINTER AND
DETERGENTS WITHIN
CABINET



INSTRUMENT WASHING



MIC WASHING



CONTAINER WASHING



WASHER DISINFECTOR WD7010

Washer disinfector with single or double sliding glass door. The latest generation of electronics with full touch screen technology incorporating visuals of safety and results. Also featuring our 35' thermal disinfection TURBO cycle. Optimisation of the running costs in versions with steam heating.

12
DIN

AUTOMATIC RACK
RECOGNITION



INTEGRATED
PRINTER



OP SHOES WASHING



CUSTOMISED WASHING



INSTRUMENT WASHING



WASHER DISINFECTOR WD7015

Washer disinfector with single or double sliding glass door. The latest generation of electronics with full touch screen technology guarantees impeccable performances. Capacity: up to 18 DIN baskets. Also featuring our 35' TURBO thermal disinfection cycle.

18
DIN

SPRAY ARM DETECTION
SYSTEM



USB PORT FOR SAFETY
STORAGE



WD7015TECO

A new thermal disinfector designed to have reduced cycle times and minimal impact on the environment, whilst maintaining a capacity of 18 DIN baskets. With its efficient heat exchanger, the WD7015TECO achieves optimal results, reducing energy consumption by 19.4% and water usage by 33% (for three phase cycles). The completion time of each full cycle is reduced by 10%, thereby increasing the CSSD's* daily washing capacity. * (Central Sterile Services Department).

ROBOTIC/DA VINCI SET
WASHING



CONTAINER WASHING



ANESTHESIA SET WASHING



SMEG S.p.A.

Founded in 1948, SMEG, whose name is the acronym of “Smalterie Metallurgiche Emiliane Guastalla”, is a proudly Italian company now in the third generation of the Bertazzoni family. Renowned for the responsible production of creations that combine technology and style, SMEG manufactures high-performance products, made with quality materials for the **home, professional and medical sectors**. SMEG **production is based in Italy**, in 5 different plants and has a strong **international presence** with direct commercial subsidiaries and with a widespread worldwide network of distributors.

SMEG S.p.A. Group is internationally recognized as one of the leaders of **Made in Italy**.



For more information, visit: www.smeg.com



Smeg S.p.A.

Via Leonardo da Vinci, 4 - 42016 Guastalla (RE) - Italy

www.smeg-instruments.com | e-mail: instruments@smeg.it