

WD2050



SINGLE PHASE WASHER DISINFECTOR, EN 15883 APPROVED

BASIC LINE - WASHER DISINFECTOR 45 CM

[Download PDF](#)

- ✓ washing up to 95 °C - antibacterial thermal disinfection
- ✓ programs: 6 default + 2 free
- ✓ microprocessor programmable electronic control
- ✓ serial port RS232 for connection to external printer or PC
- ✓ active thermodynamic drying
- ✓ high efficiency steam condenser with cold water injection
- ✓ peristaltic pump for liquid detergent
- ✓ peristaltic pump for acid liquid neutralizer
- ✓ washing chamber and internal door in stainless steel AISI 316L
- ✓ external panels in stainless steel AISI 304
- ✓ 2 independent washing levels through telescopic guides
- ✓ external dimensions LxDxH: 450 x 640 x 850 mm
- ✓ internal working dimensions LxDxH: 370 x 480 x 570 mm
- ✓ incorporated water softener with salt sensor
- ✓ cold water connection
- ✓ electronic door lock system for user safety
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 3,3 kW max
- ✓ 60 Hz version available
- ✓ Schuko plug 16 A / 250 V included
- ✓ without trolleys

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it

