


WD2145



WASHER DISINFECTOR, EN 15883 APPROVED

BASIC LINE - WASHER DISINFECTOR 45 CM

[Download PDF](#)

- ✓ washing up to 95 °C - antibacterial thermal disinfection
- ✓ programs: 9 default + 6 custom
- ✓ microprocessor programmable electronic control
- ✓ serial port RS232 for connection to external printer or PC
- ✓ active thermodynamic drying
- ✓ high efficiency steam condenser Ecoslim, with cold water injection
- ✓ peristaltic pump for liquid detergent and level sensor
- ✓ peristaltic pump for acid liquid neutralizer and level sensor
- ✓ washing chamber and internal door in stainless steel AISI 316L
- ✓ external panels in stainless steel AISI 304
- ✓ n. 2 independent washing levels through telescopic guides
- ✓ external dimensions LxDxH: 450 x 620 x 850 mm
- ✓ internal working dimensions LxDxH: 370 x 480 x 570 mm
- ✓ incorporated water softener with salt sensor
- ✓ cold water and demineralized water connections
- ✓ electronic door lock system for user safety
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 2,8 kW max
- ✓ 60 Hz version available
- ✓ without trolleys

